




The Biotechnology Education Company ®



EDVO-Kit #  
**302**

## Purification of the Restriction Enzyme *Eco* RI

**Storage: See Page 3 for  
specific storage instructions**

### **EXPERIMENT OBJECTIVE:**

In this experiment, students will purify a restriction endonuclease, test its enzyme activity, and visualize the test results by agarose gel electrophoresis.

This experiment is designed for DNA staining with InstaStain® Ethidium Bromide.

EDVOTEK, Inc. • 1-800-EDVOTEK • [www.edvotek.com](http://www.edvotek.com)

## Table of Contents

	Page
Experiment Components	3
Experiment Requirements	3
Background Information	5
Experiment Procedures	
Experiment Overview and General Instructions	7
Partial Purification of <i>Eco</i> RI	8
Analysis of <i>Eco</i> RI Activity (First Assay)	10
Quantification of <i>Eco</i> RI Activity (Second Assay)	14
Activity Determination in Units	16
Study Questions	17
Instructor's Guidelines	
Notes to the Instructor	20
Pre-Lab Preparations	23
Experiment Results and Analysis	25
Study Questions and Answers	26
Appendices	
A. 0.8% Agarose Gel Preparation Reference Tables for DNA Staining with InstaStain® Ethidium Bromide	28
B. 0.8% Agarose Gel Electrophoresis Reference Tables Quantity Gel Preparation	29
C. Staining and Visualization of DNA with InstaStain® Ethidium Bromide Cards	30
Material Safety Data Sheets	31



EVT 2010-08-17K

EDVOTEK, The Biotechnology Education Company, and InstaStain are registered trademarks of EDVOTEK, Inc.. Ready-to-Load and UltraSpec-Agarose are trademarks of EDVOTEK, Inc.

## Components & Requirements

### Experiment Components

This experiment is designed for DNA staining with InstaStain® Ethidium Bromide.

This experiment is designed for 5 laboratory groups.

<u>Components</u>	<u>Storage</u>
<b>A</b> <i>E. coli</i> RY ( <i>Eco</i> RI) Extract (lyophilized)	Refrigerator
<b>B</b> DEAE-Cellulose	Room temperature
<b>C</b> 10x Equilibration Buffer	Freezer
<b>D</b> 50% Glycerol	Freezer
<b>E</b> KCl	Room temperature
<b>F</b> <i>Eco</i> RI Reaction Buffer	Freezer
<b>G</b> Qualified Water	Freezer
<b>H</b> Lambda DNA	Freezer
<b>I</b> Lambda/ <i>Eco</i> RI Marker	Freezer
<b>J</b> <i>Eco</i> RI Dilution Buffer	Freezer

### Reagents & Supplies

- UltraSpec-Agarose™ powder
- Concentrated electrophoresis buffer
- 10x Gel Loading Solution
- Chromatography columns
- InstaStain® Ethidium Bromide

### STORAGE OF PERISHABLES

This experiment includes perishable components which were sent on wet ice. Store these components at -20°C (-4°F). Please note what type of freezer you have and store components accordingly.

#### Frost-free Freezer

Most refrigerator/freezers in homes are frost free. This means the freezer goes through warming cycles to eliminate frost (defrost cycle). If using this type of freezer, keep the enzymes in the foam chest (with the ice brick) in which they were sent. This will help maintain the enzymes at -20°C when the freezer goes through the defrost cycle.

#### Non Frost-free Freezer

These older model freezers, which are still sold but are harder to find, do not go through warming cycles. Therefore, ice will build up on freezer walls over time. If using this type of freezer, check to make sure that it maintains temperature at -20°C.

All components are intended for educational research only. They are not to be used for diagnostic or drug purposes, nor administered to or consumed by humans or animals.

THIS EXPERIMENT DOES NOT CONTAIN HUMAN DNA. None of the experiment components are derived from human sources.

### Requirements

- Horizontal gel electrophoresis apparatus
- D.C. power supply
- Automatic micropipets with tips
- Balance
- Water bath
- Ring stand and clamps
- 13 x 100 mm glass test tubes
- Assorted laboratory glassware
- Permanent markers and tape
- 1.5 ml microtest tubes
- Microwave, hot plate or burner
- Pipet pump
- 250 ml flasks or beakers
- Hot gloves
- Safety goggles and disposable laboratory gloves
- Distilled or deionized water
- UV Transilluminator

#### Technical Service Department

Mon - Fri

9:00 am to 6:00 pm ET

**I-800-EDVOTEK**

(1-800-338-6835)

FAX: (301) 340-0582

Web: [www.edvotek.com](http://www.edvotek.com)

email: [info@edvotek.com](mailto:info@edvotek.com)

EDVO-TECH SERVICE

Mon - Fri 9 am - 6 pm ET

Please have the following  
information ready:

- Experiment number and title
- Kit lot number on box or tube
- Literature version number  
(in lower right corner)
- Approximate purchase date

Online Ordering

now available

[www.edvotek.com](http://www.edvotek.com)

Visit our web site for information  
about EDVOTEK's complete line  
of "hands-on" experiments for  
biotechnology and biology education.



EVT 2010-08-17K

## Purification of the Restriction Enzyme *Eco* RI

Sequence-specific, or Type II, endonucleases are commonly known as restriction enzymes. In contrast with nonspecific endonucleases, these enzymes generate reproducible fragments from specific DNAs. They cleave double-stranded DNA by hydrolyzing two phosphodiester bonds (one per strand) within defined nucleotide sequences. Over 3,000 enzymes have been discovered since the first report by H.O. Smith and collaborators. These enzymes are extracted from a variety of bacterial strains.

Restriction Enzyme	Recognition Site
<i>Bam</i> HI <i>Bacillus amyloliquefaciens</i> H	$\begin{array}{c} \downarrow \\ 5\text{'-GGATCC-3'} \\ 3\text{'-CCTAGG-5'} \\ \uparrow \end{array}$
<i>Bgl</i> I <i>Escherichia coli</i> RY13	$\begin{array}{c} \downarrow \\ 5\text{'-GCCNNNNNGGC-3'} \\ 3\text{'-CGGNNNNNCCG-5'} \\ \uparrow \end{array}$
<i>Eco</i> RI <i>Bacillus globigii</i>	$\begin{array}{c} \downarrow \\ 5\text{'-GAATTC-3'} \\ 3\text{'-CTTAAG-5'} \\ \uparrow \end{array}$
<i>Hae</i> III <i>Haemophilus aegyptius</i>	$\begin{array}{c} \downarrow \\ 5\text{'-GGCC-3'} \\ 3\text{'-CCGG-5'} \\ \uparrow \end{array}$
<i>Hind</i> III <i>Haemophilus influenzae</i> R <sub>4</sub>	$\begin{array}{c} \downarrow \\ 5\text{'-AAGCTT-3'} \\ 3\text{'-TTCGAA-5'} \\ \uparrow \end{array}$

Figure 1: Examples of Restriction Enzymes and their recognition sites

The name of a restriction enzyme is derived from the genus and species of bacterium from which it is isolated. The first letter of the genus name and first two letters of the species are combined to form the enzyme name. This is followed by a strain designation if applicable. In many instances, a bacterial strain contains more than one restriction endonuclease. When this occurs, each enzyme is assigned a Roman numeral. For example, *Bam* HI was the first enzyme activity reported from *Bacillus amyloliquefaciens* strain H.

Most restriction enzymes are composed of two polypeptides of equal subunits with molecular weights of 20,000-25,000 or single polypeptides with molecular weights of 30,000-35,000. Enzyme activities can be differentiated from each other by their characteristic digestion patterns of small viral DNAs. The DNA from bacteriophage lambda is the most widely used substrate for screening restriction enzymes. Because it is often difficult to determine a characteristic pattern from a lambda digest, smaller DNAs, such as the replicative form of bacteriophage  $\phi$ X174 and SV40 DNA are also used as substrates. The resulting DNA restriction enzyme digests are displayed on agarose gels and visualized by staining with ethidium bromide.

A given recognition sequence in DNA can often be cleaved by more than one restriction enzyme. The term "isoschizomers" describes a group of restriction enzymes that recognize the same sequence in DNA. The sequences recognized by these enzymes are for the most part centrosymmetric "palindromic" sequences that are usually hexamers, pentamers, or tetramers. Several Type II restriction enzymes recognize DNA at a specific site and hydrolyze phosphodiester bonds at a defined distance from that site. An example of this group of enzymes is *Bgl* I, which recognizes a sequence containing two groups of specified residues separated by completely unspecified residues - GCCNNNNNGGC; it therefore generates DNA fragments with variable end groups.

## Purification of the Restriction Enzyme *Eco* RI

There is considerable diversity in the fragment termini produced in cleavage by Type II endonucleases that recognize and cleave within the same sequence. In some cases, the 5' extension may be as short as two nucleotides or as long as five. Points of cleavage on each strand may be opposite each other; this results in blunt (square ends). Several restriction endonucleases produce 3' extensions of two to four nucleotides. However, all Type II endonucleases produce fragments with a 5'-terminal phosphate and a 3'-terminal hydroxyl residue (Figure 1).

Enzymes in the Type II restriction enzyme family are amenable to purification by chromatographic procedures. Ion exchangers at nearly neutral pH are used as separation matrices after extracts have been freed of cellular nucleic acids. At this stage of purification, short-term assays often make it possible to visualize enzyme fractions that contain restriction enzymes. A variety of enzymes have been fractionated with affinity chromatography. This method takes advantage of biospecific interactions not offered by conventional fractionation methods. The advantages of affinity chromatography are speed of purification and often protection against denaturation during fractionation.

### Effects of Reaction Conditions on Restriction Enzymes

Several reports have described apparent changes in specificity of restriction endonucleases in association with altered reaction environments. Conditions that alter specificity have included changes in ionic concentration, pH of the reaction buffer, and the amounts of glycerol in the storage and the reaction mixture. For example, when lambda DNA is incubated with *Eco* RI or *Bam* HI in the presence of glycerol at various concentrations, a progressive change in the DNA digestion pattern is observed.

A change in recognition specificity of enzymes include *Bam* HI and *Eco* RI activity. The second activity is designated as ".1" (as *Bam* HI.1). A similar activity is displayed by *Eco* RI. Increasing the pH of the reaction from 7.0 to 9.0 in the absence of monovalent cations stimulates alternate activities. Decreases in the ionic strength have a similar effect.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1996, 1998, 1999, 2001, 2002, 2006, 2007, 2010 EDVOTEK, Inc., all rights reserved  
EVT 2010-08-17K

## Experiment Overview and General Instructions

### EXPERIMENT OBJECTIVE:

In this experiment, students will purify a restriction endonuclease, test its enzyme activity, and visualize the test results by agarose gel electrophoresis.



### LABORATORY SAFETY

1. Gloves and goggles should be worn routinely as good laboratory practice.
2. Exercise extreme caution when working with equipment that is used in conjunction with the heating and/or melting of reagents.
3. DO NOT MOUTH PIPET REAGENTS - USE PIPET PUMPS.
4. Exercise caution when using any electrical equipment in the laboratory.
5. Always wash hands thoroughly with soap and water after handling reagents or biological materials in the laboratory.
6. This experiment utilizes InstaStain® Ethidium Bromide for staining and visualization of DNA after gel electrophoresis. Always wear gloves when handling InstaStain® cards. Although there is only a very small amount of Ethidium bromide on InstaStain® EtBr cards, it is a listed mutagen. Wear UV-resistant safety goggles when working with ultraviolet light since it can cause irreparable damage to the eyes. Exposure to skin should also be avoided.



## Partial Purification of *Eco* RI

### Note

Loading the column and subsequent elution will be done at room temperature. The elution buffers and the fractions collected should be stored on ice as they elute from the column.

### PACKING AND EQUILIBRATING THE COLUMN

1. Vertically mount the column on a ring stand. Make sure it is straight.
2. Slide the cap onto the spout at the bottom of the column.
3. Mix the DEAE-Cellulose (ion-exchanger matrix) thoroughly by swirling or gently stirring.
4. Carefully pipet the mixed DEAE-Cellulose into the column by letting it stream down the inside walls of the column.  
  
If the flow is stopped by an air pocket, stop adding the DEAE-Cellulose and firmly tap the column until the air is removed and the exchanger flows down. Continue adding the exchanger.
5. Place an empty beaker under the column to collect wash material.
6. Remove the cap from the bottom of the column and allow the matrix to pack into the column.
7. Wash the packed column with 25 ml of Eq (1x equilibration buffer).

**Do not allow the column to dry.**



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1996, 1998, 1999, 2001, 2002, 2006, 2007, 2010 EDVOTEK, Inc., all rights reserved  
EVT 2010-08-17K



Partial Purification of *Eco* RI

## COLLECTING COLUMN FRACTIONS

1. Label eight pieces of tape 2-9 with a permanent lab marker and adhere to 13 x 100 mm test tubes. The chart to the left indicates which tubes to use for the different fractions.
2. Add 3 ml of distilled water to a test tube and use this as a reference guide for collecting the eluted fractions.
3. Slowly load the column with 1 ml of *E. coli* RY extract. Allow the extract to completely enter the column.
4. Slowly add 6 ml of Eq to the column to remove protein that is in the flow through. Collect two fractions (3 ml each) into the tubes labeled 2 and 3 and store on ice.
5. Sequentially elute the column with the following buffers. In each case, collect 3 ml fractions into the appropriate tubes and store fractions on ice immediately upon collection. Do not let the column go dry.
  - 6 ml of 0.1 M KCl. Collect two fractions, 3 ml each, into tubes 4 and 5. Store on ice.
  - 6 ml of 0.2 M KCl. Collect two 3 ml fractions into tubes 6 and 7. Store on ice.
  - 6 ml of 0.5 M KCl. Collect two 3 ml fractions into tubes 8 and 9. Store on ice.

**Table 1: Key for Identifying Fractions**

Tube	Fraction
2	(no salt)
3	(no salt)
4	0.1 M KCl
5	0.1 M KCl
6	0.2 M KCl
7	0.2 M KCl
8	0.5 M KCl
9	0.5 M KCl



## OPTIONAL STOPPING POINT

If time does not permit you to continue with *Eco* RI Activity analysis, you may freeze the fractions at  $-20^{\circ}\text{C}$  and perform the assays at a later time. Thaw the fractions at room temperature and immediately place on ice. Continue with the analysis of *Eco* RI Activity.

### Analysis of *Eco* RI Activity (First Assay)

Lambda DNA will be incubated with the fractions collected and the samples will be electrophoresed in an agarose gel to determine the peak activity of *Eco* RI endonuclease. Lambda DNA cut with *Eco* RI yields a characteristic and recognizable fragmentation pattern.

#### Important Note

Reagents listed in Table 2 will be used for two assays. Label your set of reagent tubes with your initials or group number and store reagents in the refrigerator between assays.



**Table 2: Reagents for the Incubation of *Eco* RI with Lambda DNA**

Water (G)	1 ml	on ice
<i>Eco</i> RI Rxn Buffer (F)	100 $\mu$ l	on ice
Lambda DNA (H)	100 $\mu$ l	on ice
10x Gel Load	100 $\mu$ l	
Marker	45 $\mu$ l	

1. With a permanent marker, label 9 microtest tubes 1-9. Put your initials and group number on each tube.
2. Each group will assay *Eco* RI using 2  $\mu$ l, 4  $\mu$ l, 6  $\mu$ l, 8  $\mu$ l, or 10  $\mu$ l as assigned by your instructor. In Table 3 on the next page, the "x" equals your assigned volume for analysis.
3. Use an automatic micropipet to add (40-x  $\mu$ l) of Qualified water to each of the 9 tubes
4. Use an automatic micropipet to dispense 5  $\mu$ l of the *Eco* RI Rxn Buffer and 5  $\mu$ l of Lambda DNA to each of the 9 tubes.
5. Use a clean pipet tip for each fraction and add 2  $\mu$ l, 4  $\mu$ l, 6  $\mu$ l, 8  $\mu$ l, or 10  $\mu$ l as assigned, from each fraction to the appropriate tube.

6. Cap the tubes tightly and tap on the lab bench to collect samples at the bottom of the tubes or quick spin balanced tubes in a microcentrifuge.
7. Mix the samples and incubate in a 37°C waterbath for 15 minutes.
8. Make sure your set of reagent tubes are labeled with your initials or group number and store in the refrigerator for later use in the second assay.
9. After the 15 minute incubation is complete, add 5  $\mu$ l of 10x gel loading solution to each tube to stop the reactions.

This prepares the *Eco* RI digestion products for separation by agarose gel electrophoresis.



Analysis of *Eco* RI Activity (First Assay)

## Sequence for Restriction Enzyme Reactions

Rxn Tube	Qualified Water ( $\mu$ l)	<i>Eco</i> RI Reaction Buffer ( $\mu$ l)	Lambda DNA ( $\mu$ l)	Fraction	Reaction Volume ( $\mu$ l)	37°C Incubation (minutes)	10x Gel Load ( $\mu$ l)
1	40	5	5	None	50	15	5
2	(40 - x)	5	5	x $\mu$ l tube 2 (no salt)	50	15	5
3	(40 - x)	5	5	x $\mu$ l tube 3 (no salt)	50	15	5
4	(40 - x)	5	5	x $\mu$ l tube 4 (0.1 M KCl)	50	15	5
5	(40 - x)	5	5	x $\mu$ l tube 5 (0.1 M KCl)	50	15	5
6	(40 - x)	5	5	x $\mu$ l tube 6 (0.2 M KCl)	50	15	5
7	(40 - x)	5	5	x $\mu$ l tube 7 (0.2 M KCl)	50	15	5
8	(40 - x)	5	5	x $\mu$ l tube 8 (0.5 M KCl)	50	15	5
9	(40 - x)	5	5	x $\mu$ l tube 9 (0.5 M KCl)	50	15	5

\* Volumes of *Eco* RI in fractions should be varied among different groups within the range of 2 to 5  $\mu$ l, with 1  $\mu$ l increments. Water in the assay should be adjusted accordingly.

\*\* To be added after incubation at 37°C.



## OPTIONAL STOPPING POINT

If time does not permit you to continue with agarose gel electrophoresis at this time, you may freeze the fractions at -20°C and perform the electrophoresis at a later date. Thaw the fractions at room temperature and heat the samples at 65°C before loading the gel.

## Analysis of *Eco* RI Activity (First Assay)

If you are unfamiliar with agarose gel preparation and electrophoresis, detailed instructions and helpful resources are available at [www.edvotek.com](http://www.edvotek.com)

### Important Note



Continue heating until the final solution appears clear (like water) without any undissolved particles. Check the solution carefully. If you see "crystal" particles, the agarose is not completely dissolved.

### AGAROSE GEL REQUIREMENTS FOR THE FIRST ASSAY

- Recommended gel size: 7 x 14 cm
- Number of sample wells required: 10
- Agarose gel concentration: 0.8%

### PREPARING THE AGAROSE GEL

1. Close off the open ends of a clean and dry gel bed (casting tray) by using rubber dams or tape.
2. Place a well-former template (comb) in the first set of notches at the end of the bed. Make sure the comb sits firmly and evenly across the bed.
3. To a 250 ml flask or beaker, add agarose powder and buffer as indicated in the Reference Tables (Appendix A) provided by your instructor. Swirl the mixture to disperse clumps of agarose powder.
4. With a marking pen, indicate the level of the solution volume on the outside of the flask.
5. Heat the mixture using a microwave oven or burner to dissolve the agarose powder.
6. Cool the agarose solution to 60°C with careful swirling to promote even dissipation of heat. If detectable evaporation has occurred, add distilled water to bring the solution up to the original volume marked in step 4.

### After the gel is cooled to 60°C:

7. Place the bed on a level surface and pour the cooled agarose solution into the bed.
8. Allow the gel to completely solidify. It will become firm and cool to the touch after approximately 20 minutes.
9. After the gel is solidified, be careful not to damage or tear the wells while removing the rubber dams or tape and comb(s) from the gel bed.
10. Place the gel (on its bed) into the electrophoresis chamber, properly oriented, centered and level on the platform.
11. Fill the electrophoresis apparatus chamber with the appropriate amount of diluted (1x) electrophoresis buffer (refer to Table B on the instruction sheet from the Appendix provided by your instructor).



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1996, 1998, 1999, 2001, 2002, 2006, 2007, 2010 EDVOTEK, Inc., all rights reserved  
EVT 2010-08-17K

## Analysis of *Eco* RI Activity (First Assay)



### LOADING THE SAMPLES

This experiment is designed for staining with InstaStain® Ethidium Bromide. The amount of sample that should be loaded is 18-20  $\mu$ l. Make sure the gel is completely submerged under buffer before loading the samples and conducting electrophoresis.

Lane	Tube	
1	Marker	Lambda Eco RI Marker
2	1	Uncut Lambda DNA
3	2	Lambda + 2 (no salt)
4	3	Lambda + 3 (no salt)
5	4	Lambda + 4 (0.1 M KCl)
6	5	Lambda + 5 (0.1 M KCl)
7	6	Lambda + 6 (0.2 M KCl)
8	7	Lambda + 7 (0.2 M KCl)
9	8	Lambda + 8 (0.5 M KCl)
10	9	Lambda + 9 (0.5 M KCl)

#### Reminder:

Before loading the samples, make sure the gel is properly oriented in the apparatus chamber.



Electrophoresis can be completed in 15-20 minutes under optimal conditions. For Time and Voltage recommendations, refer to Table C (from Appendix A or B).

### RUNNING THE GEL

1. After the DNA samples are loaded, properly orient the cover and carefully snap it onto the electrode terminals.
2. Insert the plugs of the black and red wires into the corresponding inputs of the power source.
3. Set the power source at the required voltage and conduct electrophoresis for the length of time determined by your instructor.
4. Check to see that current is flowing properly - you should see bubbles forming on the two platinum electrodes.
5. After the electrophoresis is completed, disconnect the power and remove the gel from the bed for staining.

### STAINING AND VISUALIZATION OF DNA

After electrophoresis, agarose gels require staining to visualize the separated DNA samples. Your instructor will provide instructions for DNA staining with InstaStain® Ethidium Bromide.

### Quantification of *Eco* RI Activity (Second Assay)

Units of enzyme activity are defined by convention. A restriction enzyme unit is defined as the amount of enzyme activity that will digest 1  $\mu\text{g}$  of lambda DNA at 37°C within one hour under the defined assay conditions. To determine the total units of *Eco* RI purified in this experiment, re-assay pooled enzyme fractions at various enzyme dilutions to determine the minimum amount of enzyme that yields complete digestion of 1  $\mu\text{g}$  of lambda DNA.

**Table 4: Dilution of Pooled *Eco* RI**

Pooled Enzyme ( $\mu\text{l}$ )	Dilution Buffer ( $\mu\text{l}$ )	Total Volume ( $\mu\text{l}$ )	Dilution Factor
10	0	10	0
5	5	10	1:2
10	20	30	1:3
10	30	40	1:4
10	90	100	1:10
5	95	100	1:20

1. Pool the enzyme fractions that have *Eco* RI activity as judged by the first assay. If a fraction has only a trace of activity, do not pool it since it will dilute the enzyme which may result in activity loss.
2. Measure and record the volume of pooled *Eco* RI fractions.
3. Gently mix the pooled fraction to get a representative sample for assaying.
4. Dilute the pooled *Eco* RI enzyme fraction with *Eco* RI Dilution Buffer, using the dilution factors indicated in Table 4.
5. Prepare each *Eco* RI dilution (from Table 4) for incubation as outlined in Table 5. Use the remaining reagents from the first assay that were stored in the refrigerator.
6. After completing the incubations as outlined in Table 5, add 10x Gel Loading Solution to each tube to stop the reactions.
7. Separate the *Eco* RI digestion products by agarose gel electrophoresis.

Store all fractions on ice.

Label fractions according to dilution factors. 10  $\mu\text{l}$  of each dilution will be used as shown in Table 5.

**Table 5: Assay to Determine Total Units of *Eco* RI**

Rxn Tube	Qualified Water ( $\mu\text{l}$ )	<i>Eco</i> RI Reaction Buffer ( $\mu\text{l}$ )	Lambda DNA ( $\mu\text{l}$ )	<i>Eco</i> RI Dilution (from Table 4)	Reaction Volume ( $\mu\text{l}$ )	37°C Incubation (minutes)	10x Gel Load* ( $\mu\text{l}$ )
1	40	5	5	None	50	30	5
2	30	5	5	10 $\mu\text{l}$ of 0	50	30	5
3	30	5	5	10 $\mu\text{l}$ of 1:2	50	30	5
4	30	5	5	10 $\mu\text{l}$ of 1:3	50	30	5
5	30	5	5	10 $\mu\text{l}$ of 1:4	50	30	5
6	30	5	5	10 $\mu\text{l}$ of 1:10	50	30	5
7	30	5	5	10 $\mu\text{l}$ of 1:20	50	30	5

\*To be added after 37°C incubation



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1996, 1998, 1999, 2001, 2002, 2006, 2007, 2010 EDVOTEK, Inc., all rights reserved  
EVT 2010-08-17K

## Quantification of *Eco* RI Activity (Second Assay)

### AGAROSE GEL REQUIREMENTS FOR THE SECOND ASSAY

- Recommended gel size: 7 x 14 cm
- Number of sample wells required: 8
- Agarose gel concentration: 0.8%

1. Prepare a 0.8% agarose gel for the second assay according to instructions previously described.
2. Load 20  $\mu$ l of each DNA sample in the following manner:

Lane	Tube	
1	Marker	Lambda <i>Eco</i> RI Marker
2	1	Uncut Lambda DNA
3	2	Lambda + Undiluted <i>Eco</i> RI
4	3	Lambda + 1:2 Dilution
5	4	Lambda + 1:3 Dilution
6	5	Lambda + 1:4 Dilution
7	6	Lambda + 1:10 Dilution
8	7	Lambda + 1:20 Dilution

3. After the samples are loaded, conduct electrophoresis and stain the gel with InstaStain® Ethidium Bromide for visualization.
3. Examine the gel or take a photograph to determine which lane gives complete digestion determined as follows:
  - No undigested or partially digested lambda DNA is visible.
  - All the DNA digestion products (5 bands) are visible.

## Activity Determination in Units

Restriction enzyme unit = amount of enzyme activity that will digest 1  $\mu$ g of lambda DNA at 37°C within one hour.

### DETERMINATION OF TOTAL ACTIVITY

Total units (units) is the amount of enzyme activity recovered from the preparation. It does not indicate the level of enzyme purity.

$$\text{Total Activity (units)} = \frac{\text{Pooled volume } (\mu\text{l})}{\text{Volume used for assay } (\mu\text{l})} \times \text{Dilution factor}$$

#### Example for Determining Total units.

Pooled volume is 9 ml = 9000  $\mu$ l  
*Eco* RI volume for assay = 10  $\mu$ l  
 Dilution factor = 4

$$\frac{9000 \mu\text{l}}{10 \mu\text{l}} \times 4 = 3600 \text{ units}$$

for a 30 minute digestion

Conversion for a 1 hour digestion assay:

$$\text{Total Activity units} = 3600 \text{ units} \times 2 = 7200$$

### SPECIFIC ACTIVITY DETERMINATION (OPTIONAL)

Specific activity is defined as the number of enzyme units per mg of total protein in the enzyme fraction. The less total protein the *Eco* RI fraction contains, the higher is its specific activity.

- For this experiment we have equated 1.0 absorbance unit at A<sub>280</sub>. In 9 ml, the amount of protein is 0.2 mg/ml  $\times$  9 ml = 1.8 ml.

$$\text{Specific Activity} = \frac{\text{Total units}}{\text{mg of protein}}$$

#### Example for Determining Specific Activity

Total units: 7200 units for the total volume of 9 ml  
 Total mg. of protein - 1.8 mg

$$\text{Specific Activity} = \frac{7200 \text{ units}}{1.8 \text{ mg}} = 4,000 \text{ units/mg}$$



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1996, 1998, 1999, 2001, 2002, 2006, 2007, 2010 EDVOTEK, Inc., all rights reserved  
 EVT 2010-08-17K



## Study Questions

Answer the following study questions in your laboratory notebook or on a separate worksheet.

1. What is the recognition site for *Eco* RI?
2. How is *E. coli* host DNA protected against action of the *Eco* RI endonuclease?
3. How many *Eco* RI sites are there in lambda DNA?
4. What is the difference between total activity versus specific activity?



EVT 2010-08-17K

## Instructor's Guide

Class size, length of laboratory sessions, and availability of equipment are factors which must be considered in the planning and the implementation of this experiment with your students. These guidelines can be adapted to fit your specific set of circumstances. If you do not find the answers to your questions in this section, a variety of resources are continuously being added to the EDVOTEK web site. In addition, Technical Service is available from 9:00 am to 6:00 pm, Eastern time zone. Call for help from our knowledgeable technical staff at 1-800-EDVOTEK (1-800-338-6835).

### EDUCATIONAL RESOURCES

#### Electrophoresis Hints, Help and Frequently Asked Questions

EDVOTEK Experiments are designed for maximum success in the classroom setting. However, even the most experienced students and teachers occasionally encounter experimental problems or difficulties. The EDVOTEK web site provides several suggestions and reminders for conducting electrophoresis, as well as answers to frequently asked electrophoresis questions.



**EDVO-TECH SERVICE**  
**1-800-EDVOTEK**  
(1-800-338-6835)  
Mon - Fri 9 am - 6 pm ET

**Technical Service Department**  
Mon - Fri  
9:00 am to 6:00 pm ET  
FAX: (301) 340-0582  
Web: [www.edvotek.com](http://www.edvotek.com)  
email: [info@edvotek.com](mailto:info@edvotek.com)

Please have the following information ready:

- Experiment number and title
- Kit lot number on box or tube
- Literature version number (in lower right corner)
- Approximate purchase date

Online Ordering  
now available



Visit our web site for information about EDVOTEK's complete line of "hands-on" experiments for biotechnology and biology education.

## Notes to the Instructor:

**MICROPIPETTING BASICS AND PRACTICE GEL LOADING**

Accurate pipeting is critical for maximizing successful experiment results. EDVOTEK Series 300 experiments are designed for students who have had previous experience with agarose gel electrophoresis and micropipetting techniques. If your students are unfamiliar with using micropipets, EDVOTEK highly recommends that students perform Experiment # S-44, Micropipetting Basics, or other Series 100 or 200 electrophoresis experiment prior to conducting this advanced level experiment.

**APPROXIMATE TIME REQUIREMENTS**

- **Pre-lab preparations**  
Pre-lab preparations and dispensing of biologicals and reagents take approximately 1-2 hours.
- **Restriction Enzyme Digestion**  
The approximate time required for students to perform the restriction enzyme digestion and prepare samples for electrophoresis is 50-75 minutes. Extending the restriction enzyme digest incubation time to 60 minutes will help ensure complete cleavage of DNA.
- **Agarose Gel preparation**  
Whether you choose to prepare the gel(s) in advance or have the students prepare their own, allow approximately 30-40 minutes for this procedure. Generally, 20 minutes of this time is required for gel solidification. See section "Options for Preparing Agarose Gels" below.

Volts	EDVOTEK Electrophoresis Model	
	M6+ Minimum / Maximum	M12 & M36 Minimum / Maximum
150	15 / 20 min	25 / 35 min
125	20 / 30 min	35 / 45 min
70	35 / 45 min	60 / 90 min
50	50 / 80 min	95 / 130 min

- **Conducting Electrophoresis**  
The approximate time for electrophoresis will vary from 15 minutes to 2 hours. Generally, the higher the voltage applied, the faster the samples migrate. However, depending upon the apparatus configuration and the distance between the two electrodes, individual electrophoresis units will separate DNA at different rates. Follow manufacturer's recommendations. Time and Voltage recommendations for EDVOTEK equipment are outlined in Table C.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1996, 1998, 1999, 2001, 2002, 2006, 2007, 2010 EDVOTEK, Inc., all rights reserved  
EVT 2010-08-17K

**Notes to the Instructor:****OPTIONS FOR PREPARING AGAROSE GELS**

This experiment is designed for DNA staining after electrophoresis with InstaStain® Ethidium Bromide. There are several options for preparing agarose gels for the experiment.

1. Individual Gel Casting:  
Each student lab group can be responsible for casting their own individual gel prior to conducting the experiment.
2. Preparing Gels in Advance:  
Gels may be prepared ahead and stored for later use. Solidified gels can be stored under buffer in the refrigerator for up to 2 weeks.

Do not store gels at -20°C. Freezing will destroy the gels.

Gels that have been removed from their trays for storage, should be "anchored" back to the tray with a few drops of hot, molten agarose before placing the gels into the apparatus for electrophoresis. This will prevent the gels from sliding around in the trays and the chambers.

3. Batch Gel Preparation:  
A batch of agarose gel can be prepared for sharing by the class. To save time, a larger quantity of UltraSpec-Agarose can be prepared for sharing by the class. See instructions for "Batch Gel Preparation".

**GEL CONCENTRATION AND VOLUME**

The gel concentration required is 0.8%. Prepare 7 x 14 cm gels according to Table A.1 or A.2 in Appendix A.

**Notes to the Instructor:****GEL STAINING AND DESTAINING AFTER ELECTROPHORESIS**

This experiment features InstaStain® Ethidium Bromide for gel staining after electrophoresis. It is a proprietary staining method which saves time and reduces liquid waste. DNA staining with InstaStain® Methylene Blue is not recommended because it will not yield optimal results.

**Instastain® Ethidium Bromide**

- InstaStain® Ethidium Bromide                      Appendix C

Optimal visualization of DNA fragments on gels is obtained by staining with InstaStain® Ethidium Bromide (InstaStain® EtBr) cards.

Caution: Ethidium Bromide is a listed mutagen. Disposal of the InstaStain® EtBr cards, which contain only a few micrograms of ethidium bromide, is minimal compared to the large volume of liquid waste generated by traditional ethidium bromide staining procedures. Disposal of InstaStain® cards and gels should follow institutional guidelines for chemical waste.

**PHOTODOCUMENTATION OF DNA (OPTIONAL)**

There are many different photodocumentation systems available, including digital systems that are interfaced directly with computers. Specific instructions will vary depending upon the type of photodocumentation system you are using.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1996, 1998, 1999, 2001, 2002, 2006, 2007, 2010 EDVOTEK, Inc., all rights reserved  
EVT 2010-08-17K

## Pre-Lab Preparations

PARTIAL PURIFICATION OF *ECO* RI  
(Packing the Column and Collecting Fractions)

The 10x equilibration buffer used to hydrate the DEAE-Cellulose contains potassium phosphate, pH 7.4, EDTA, and  $\beta$ -mercaptoethanol.

## DEAE-Cellulose Matrix

- Hydrate the ion-exchanger, DEAE-Cellulose (B) in 35 ml of 10x equilibration buffer (C).
- Stir occasionally for a minimum of 30 minutes.
- Aliquot 6 ml for each of the five groups.

## Quick Reference

## Summary of Reagent Preparations

DEAE-Cellulose (Matrix, B)	6 ml	on ice
Eq Buffer (1x equil buffer, diluted C)	35 ml	on ice
0.1 M KCl	6 ml	on ice
0.2 M KCl	6 ml	on ice
0.5 M KCl	6 ml	on ice
<i>E. coli</i> RY extract (A)	1 ml	on ice

## Buffers

- Prepare 500 ml of 1x equilibration buffer (Eq) in a 600 ml flask or beaker. To prepare, add the following and stir thoroughly:

350 ml	Distilled water
50 ml	10x Equilibration buffer (C)
100 ml	50% glycerol (D)

Use this prepared Eq buffer to make buffers in step 5.

## Each group requires:

- DEAE-Cellulose
- Eq Buffer
- 0.1 M KCl
- 0.2 M KCl
- 0.5 M KCl
- E. coli* RY extract
- 1 chromatography column
- 1 ring stand with clamp
- 9 test tubes (13 x 100 mm)
- 10 microtest tubes
- Automatic micropipet & tips
- 5 ml pipets and pipet pumps

- To make Eq + KCl Buffers, mix the following:

	Eq Buffer (1x)	KCl
0.1 M KCl	100 ml	0.75 g
0.2 M KCl	100 ml	1.5 g
0.5 M KCl	100 ml	3.75 g

*E. coli* Cell Extract Containing *Eco* RI Restriction Enzyme

- Re-hydrate the sample by adding 0.5 ml of distilled or deionized water to tube component A and let sit for 5 minutes.
- Mix vigorously by vortexing and transfer the entire contents to a 50 ml conical tube. Rinse tube A six times - each time with 1 ml of 1x Equilibration buffer (diluted component C) and add the rinse material to the 50 ml conical tube. Mix the tube well.
- Label 5 tubes "*E. coli* RY extract". Aliquot 1 ml of the re-hydrated extract for each of the student groups. Store extracts on ice.

## Pre-Lab Preparations

**ANALYSIS AND QUANTIFICATION OF *ECO* RI ACTIVITY  
(First and Second Assays)****Incubation of Fractions with Lambda DNA**

Important: Students should be reminded that the reagents they receive are for two assays.

1. Label 5 tubes "water" and dispense 1 ml Qualified Water (G) into the tubes. Store on ice.
2. Label 5 tubes "*Eco* RI Rxn Buffer" and dispense 100  $\mu$ l of *Eco* RI Reaction Buffer (F) into the tubes. Store on ice.

Quick Reference

**Summary of Reagent Preparations****Reagents for First & Second Assays**

Water (G)	1 ml	on ice
<i>Eco</i> RI Rxn Buffer (F)	100 $\mu$ l	on ice
Lambda DNA (H)	100 $\mu$ l	on ice
10x Gel Load	100 $\mu$ l	
Marker	45 $\mu$ l	

**Additional Reagent for Second Assay**

Lambda DNA (H)	250 $\mu$ l	on ice
----------------	-------------	--------

3. Label 5 tubes "Lambda DNA" and dispense 100  $\mu$ l of Lambda DNA (H) into the tubes. Store on ice.
4. Label 5 tubes "10x Gel Load" and dispense 100  $\mu$ l 10x Gel Loading Solution into the tubes.
5. Label 5 tubes "Marker" and dispense 45  $\mu$ l Lambda/*Eco* RI Marker (I) into the tubes.
6. Label 5 tubes "*Eco* RI Diln Buffer" and dispense 250  $\mu$ l of *Eco* RI Dilution Buffer (J) into the tubes. Store on ice.
7. Have a 37°C waterbath ready for *Eco* RI activity analysis.




Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1996, 1998, 1999, 2001, 2002, 2006, 2007, 2010 EDVOTEK, Inc., all rights reserved  
EVT 2010-08-17K



Experiment Results and Analysis

PARTIAL PURIFICATION OF *ECO* RI (FIRST ASSAY)



In the idealized schematic, the relative positions of DNA fragments are shown but are not depicted to scale. The schematic depicts an idealized gel result for identifying column fractions with *Eco* RI activity.

Lane	Tube	
1	Marker	Lambda <i>Eco</i> RI Markers
2	1	Uncut Lambda DNA
3	2	Lambda + 2 (no salt)
4	3	Lambda + 3 (no salt)
5	4	Lambda + 4 (0.1 M KCl)
6	5	Lambda + 5 (0.1 M KCl)
7	6	Lambda + 6 (0.2 M KCl)
8	7	Lambda + 7 (0.2 M KCl)
9	8	Lambda + 8 (0.5 M KCl)
10	9	Lambda + 9 (0.5 M KCl)

\* Results may vary between different groups and from the schematic depicted to the left. Some bands may be faint and thus difficult to see. You may also see extra bands due to partial digestion of the DNA. The amount of activity in the flow through may also vary.

Results of the second assay will show varying results depending upon the amount of purified enzyme activity.

**Please refer to the kit  
insert for the Answers to  
Study Questions**

## Appendices

### Appendices

- A 0.8% Agarose Gel Preparation For DNA Staining with InstaStain® Ethidium Bromide
- B 0.8% Quantity Preparations for Agarose Gel Electrophoresis
- C Staining and Visualization of DNA with InstaStain® Ethidium Bromide Cards

### Material Safety Data Sheets

### 0.8% Agarose Gel Preparation Reference Tables for DNA Staining with InstaStain® Ethidium Bromide

↓  
If preparing the gel with concentrated (50x) buffer, use Table A.1.

Table A.1 Individual 0.8%\* UltraSpec-Agarose™ Gel DNA Staining with InstaStain® EtBr

Size of Gel (cm)	Amt of Agarose + (g)	Concentrated Buffer (50X) (ml)	Distilled Water (ml)	Total Volume (ml)
7 × 7	0.2	0.5	24.5	25
7 × 14	0.4	1.0	49.0	50

\* 0.77 UltraSpec-Agarose™ gel percentage rounded up to 0.8%

↓  
If preparing the gel with diluted (1x) buffer, use Table A.2.

Table A.2 Individual 0.8%\* UltraSpec-Agarose™ Gel DNA Staining with InstaStain® Ethidium Bromide

Size of Gel (cm)	Amt of Agarose + (g)	Diluted Buffer (1x) (ml)
7 × 7	0.2	25
7 × 14	0.4	50

For DNA analysis, the recommended electrophoresis buffer is Tris-acetate-EDTA, pH 7.8. The formula for diluting EDVOTEK (50x) concentrated buffer is one volume of buffer concentrate to every 49 volumes of distilled or deionized water. Prepare buffer as required for your electrophoresis unit.

Table B Electrophoresis (Chamber) Buffer

EDVOTEK Model #	Total Volume Required (ml)	Dilution	
		50x Conc. Buffer (ml)	+ Distilled Water (ml)
M6+	300	6	294
M12	400	8	392
M36 (blue)	500	10	490
M36 (clear)	1000	20	980

Table C Time and Voltage Recommendations

Volts	EDVOTEK Electrophoresis Model	
	M6+ Minimum / Maximum	M12 & M36 Minimum / Maximum
150	15 / 20 min	25 / 35 min
125	20 / 30 min	35 / 45 min
70	35 / 45 min	60 / 90 min
50	50 / 80 min	95 / 130 min

Time and Voltage recommendations for EDVOTEK equipment are outlined in Table C. The approximate time for electrophoresis will vary from approximately 15 minutes to 2 hours depending upon various factors. Conduct electrophoresis for the length of time determined by your instructor.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1996, 1998, 1999, 2001, 2002, 2006, 2007, 2010 EDVOTEK, Inc., all rights reserved

EVT 2010-08-17K

## 0.8% Agarose Gel Electrophoresis Reference Tables Quantity Preparations

To save time, electrophoresis buffer and agarose gel solution can be prepared in larger quantities for sharing by the class. Unused diluted buffer can be used at a later time and solidified agarose gel can be remelted.

Table  
D

### Bulk Preparation of Electrophoresis Buffer

Concentrated Buffer (50x) (ml)	+	Distilled Water (ml)	=	Total Volume (ml)
60		2,940		3000 (3 L)

### BULK ELECTROPHORESIS BUFFER

Quantity (bulk) preparation for 3 liters of 1x electrophoresis buffer is outlined in Table D.

### BATCH AGAROSE GELS (0.8%)

For quantity (batch) preparation of 0.8% agarose gels, see Table E.

Table  
E


### Batch Preparation of 0.8%\* UltraSpec-Agarose™

Amt of Agarose (g)	+	Concentrated Buffer (50x) (ml)	+	Distilled Water (ml)	=	Total Volume (ml)
3.0		7.5		382.5		390

\*0.77% UltraSpec-Agarose™ gel percentage rounded up to 0.8%

Note: The UltraSpec-Agarose™ kit component is often labeled with the amount it contains. In many cases, the entire contents of the bottle is 3.0 grams. Please read the label carefully. If the amount of agarose is not specified or if the bottle's plastic seal has been broken, weigh the agarose to ensure you are using the correct amount.

1. Use a 500 ml flask to prepare the diluted gel buffer
2. Pour 3.0 grams of UltraSpec-Agarose™ into the prepared buffer. Swirl to disperse clumps.
3. With a marking pen, indicate the level of solution volume on the outside of the flask.
4. Heat the agarose solution as outlined previously for individual gel preparation. The heating time will require adjustment due to the larger total volume of gel buffer solution.
5. Cool the agarose solution to 60°C with swirling to promote even dissipation of heat. If evaporation has occurred, add distilled water to bring the solution up to the original volume as marked on the flask in step 3.
6. Dispense the required volume of cooled agarose solution for casting each gel. The volume required is dependent upon the size of the gel bed.
7. Allow the gel to completely solidify. It will become firm and cool to the touch after approximately 20 minutes. Then proceed with preparing the gel for electrophoresis.



60°C

Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1996, 1998, 1999, 2001, 2002, 2006, 2007, 2010 EDVOTEK, Inc., all rights reserved

EVT 2010-08-17K

## Staining and Visualization of DNA

### INSTASTAIN® ETHIDIUM BROMIDE CARDS

1. After electrophoresis, place the gel on a piece of plastic wrap on a flat surface. Moisten the gel with a few drops of electrophoresis buffer.
2. Wearing gloves, remove the clear plastic protective sheet, and place the unprinted side of the InstaStain® EtBr card on the gel.
3. Firmly run your fingers over the entire surface of the InstaStain® EtBr. Do this several times.



4. Place the gel casting tray and a small empty beaker on top to ensure that the InstaStain® card maintains direct contact with the gel surface.

Allow the InstaStain® EtBr card to stain the gel for 10-15 minutes.

5. After 10-15 minutes, remove the InstaStain® EtBr card. Transfer the gel to a ultraviolet (300 nm) transilluminator for viewing. Be sure to wear UV protective goggles.

Visit our web site for an animated demonstration of InstaStain® EtBr.

[www.edvotek.com](http://www.edvotek.com)

**Caution:** Ethidium Bromide is a listed mutagen.


### Disposal of InstaStain

Disposal of InstaStain® cards and gels should follow institutional guidelines for chemical waste.

### Additional Notes About Staining

- If bands appear faint, or if you are not using EDVOTEK UltraSpec-Agarose™, gels may take longer to stain with InstaStain® EtBr. Repeat staining and increase the staining time an additional 10-15 minutes.
- Gels stained alternatively with InstaStain Methylene Blue or liquid methylene blue may fade with time. Re-stain the gel to visualize the DNA bands.
- DNA 200 bp markers should be visible after staining even if the amplified DNA samples are faint or absent. If markers are not visible, troubleshoot for problems with the electrophoretic separation.

Do not stain gel(s) in the electrophoresis apparatus.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1996, 1998, 1999, 2001, 2002, 2006, 2007, 2010 EDVOTEK, Inc., all rights reserved  
EVT 2010-08-17K

The Biotechnology Education Company® • 1-800-EDVOTEK • [www.edvotek.com](http://www.edvotek.com)

**Material Safety Data Sheets**  
Full size (8.5 x 11") pdf copy of MSDS available at [www.edvotek.com](http://www.edvotek.com) or by request.

IDENTITY (As Used on Label and List)		50x Electrophoresis Buffer																	
<p><b>Material Safety Data Sheet</b> May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.</p>																			
<p><b>Section I</b> Manufacturer's Name <b>EDVOTEK, Inc.</b> Address (Number, Street, City, State, Zip Code) <b>14676 Rothgeb Drive Rockville, MD 20850</b></p>		<p>Emergency Telephone Number <b>(301) 251-5990</b> Telephone Number for information <b>(301) 251-5990</b> Date Prepared <b>10/05/06</b> Signature of Preparer (optional)</p>																	
<p><b>Section II - Hazardous Ingredients/Identify Information</b> Hazardous Components (Specific Chemical Identity, Common Name(s)) OSHA PEL ACGIH TLV Recommended % (Optional) This product contains no hazardous materials as defined by the OSHA Hazard Communication Standard.</p>																			
<p><b>Section III - Physical/Chemical Characteristics</b></p> <table border="1"> <tr> <td>Boiling Point</td> <td>No data</td> <td>Specific Gravity (at 0 = 1)</td> <td>No data</td> </tr> <tr> <td>Vapor Pressure (mm Hg)</td> <td>No data</td> <td>Melting Point</td> <td>No data</td> </tr> <tr> <td>Vapor Density (AIR = 1)</td> <td>No data</td> <td>Evaporation Rate (Butyl Acetate = 1)</td> <td>No data</td> </tr> <tr> <td>Solubility in Water</td> <td>Insoluble - cold</td> <td>Appearance and Odor</td> <td>White powder, no odor</td> </tr> </table>				Boiling Point	No data	Specific Gravity (at 0 = 1)	No data	Vapor Pressure (mm Hg)	No data	Melting Point	No data	Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data	Solubility in Water	Insoluble - cold	Appearance and Odor	White powder, no odor
Boiling Point	No data	Specific Gravity (at 0 = 1)	No data																
Vapor Pressure (mm Hg)	No data	Melting Point	No data																
Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data																
Solubility in Water	Insoluble - cold	Appearance and Odor	White powder, no odor																
<p><b>Section IV - Physical/Chemical Characteristics</b> Flash Point (Method Used) N.D. = No data Flammable Limits LEL N.D. UEL N.D.</p> <p>Extinguishing Media Water spray, dry chemical, carbon dioxide, halon or standard foam</p> <p>Special Fire Fighting Procedures Possible fire hazard when exposed to heat or flame</p> <p>Unusual Fire and Explosion Hazards None</p>																			

IDENTITY (As Used on Label and List)		50x Electrophoresis Buffer																	
<p><b>Material Safety Data Sheet</b> May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.</p>																			
<p><b>Section I</b> Manufacturer's Name <b>EDVOTEK, Inc.</b> Address (Number, Street, City, State, Zip Code) <b>14676 Rothgeb Drive Rockville, MD 20850</b></p>		<p>Emergency Telephone Number <b>(301) 251-5990</b> Telephone Number for information <b>(301) 251-5990</b> Date Prepared <b>10/05/06</b> Signature of Preparer (optional)</p>																	
<p><b>Section II - Hazardous Ingredients/Identify Information</b> Hazardous Components (Specific Chemical Identity, Common Name(s)) OSHA PEL ACGIH TLV Recommended % (Optional) This product contains no hazardous materials as defined by the OSHA Hazard Communication Standard.</p>																			
<p><b>Section III - Physical/Chemical Characteristics</b></p> <table border="1"> <tr> <td>Boiling Point</td> <td>No data</td> <td>Specific Gravity (at 0 = 1)</td> <td>No data</td> </tr> <tr> <td>Vapor Pressure (mm Hg)</td> <td>No data</td> <td>Melting Point</td> <td>No data</td> </tr> <tr> <td>Vapor Density (AIR = 1)</td> <td>No data</td> <td>Evaporation Rate (Butyl Acetate = 1)</td> <td>No data</td> </tr> <tr> <td>Solubility in Water</td> <td>Appreciable, (greater than 10%)</td> <td>Appearance and Odor</td> <td>Clear, liquid, slight vinegar odor</td> </tr> </table>				Boiling Point	No data	Specific Gravity (at 0 = 1)	No data	Vapor Pressure (mm Hg)	No data	Melting Point	No data	Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data	Solubility in Water	Appreciable, (greater than 10%)	Appearance and Odor	Clear, liquid, slight vinegar odor
Boiling Point	No data	Specific Gravity (at 0 = 1)	No data																
Vapor Pressure (mm Hg)	No data	Melting Point	No data																
Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data																
Solubility in Water	Appreciable, (greater than 10%)	Appearance and Odor	Clear, liquid, slight vinegar odor																
<p><b>Section IV - Physical/Chemical Characteristics</b> Flash Point (Method Used) N.D. = No data Flammable Limits LEL N.D. UEL N.D.</p> <p>Extinguishing Media Use extinguishing media appropriate for surrounding fire.</p> <p>Special Fire Fighting Procedures Wear protective equipment and SCBA with full facepiece operated in positive pressure mode.</p> <p>Unusual Fire and Explosion Hazards None identified</p>																			

IDENTITY (As Used on Label and List)		Enzyme Reaction Dilution Buffer																	
<p><b>Material Safety Data Sheet</b> May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.</p>																			
<p><b>Section I</b> Manufacturer's Name <b>EDVOTEK, Inc.</b> Address (Number, Street, City, State, Zip Code) <b>14676 Rothgeb Drive Rockville, MD 20850</b></p>		<p>Emergency Telephone Number <b>(301) 251-5990</b> Telephone Number for information <b>(301) 251-5990</b> Date Prepared <b>11/07/06</b> Signature of Preparer (optional)</p>																	
<p><b>Section II - Hazardous Ingredients/Identify Information</b> Hazardous Components (Specific Chemical Identity, Common Name(s)) OSHA PEL ACGIH TLV Recommended % (Optional) This product contains no hazardous materials as defined by the OSHA Hazard Communication Standard.</p>																			
<p><b>Section III - Physical/Chemical Characteristics</b></p> <table border="1"> <tr> <td>Boiling Point</td> <td>No data</td> <td>Specific Gravity (at 0 = 1)</td> <td>No data</td> </tr> <tr> <td>Vapor Pressure (mm Hg)</td> <td>No data</td> <td>Melting Point</td> <td>N/A</td> </tr> <tr> <td>Vapor Density (AIR = 1)</td> <td>No data</td> <td>Evaporation Rate (Butyl Acetate = 1)</td> <td>No data</td> </tr> <tr> <td>Solubility in Water</td> <td>soluble</td> <td>Appearance and Odor</td> <td>Clear liquid no odor, dry chemical, carbon dioxide, water, spray or foam.</td> </tr> </table>				Boiling Point	No data	Specific Gravity (at 0 = 1)	No data	Vapor Pressure (mm Hg)	No data	Melting Point	N/A	Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data	Solubility in Water	soluble	Appearance and Odor	Clear liquid no odor, dry chemical, carbon dioxide, water, spray or foam.
Boiling Point	No data	Specific Gravity (at 0 = 1)	No data																
Vapor Pressure (mm Hg)	No data	Melting Point	N/A																
Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data																
Solubility in Water	soluble	Appearance and Odor	Clear liquid no odor, dry chemical, carbon dioxide, water, spray or foam.																
<p><b>Section IV - Physical/Chemical Characteristics</b> Flash Point (Method Used) No data Flammable Limits LEL No data UEL No data</p> <p>Extinguishing Media Use extinguishing media appropriate to surrounding fire</p> <p>Special Fire Fighting Procedures Remove container from fire if possible.</p> <p>Unusual Fire and Explosion Hazards May produce toxic gases</p>																			

Section V - Reactivity Data		Section V - Reactivity Data	
Stability	Unstable X Stable	Conditions to Avoid	None
Incompatibility	No data available	Hazardous Decomposition or Byproducts	None
Hazardous Decomposition or Byproducts	None	Hazardous Polymerization	None
Hazardous Polymerization	None	Routes of Entry:	Inhalation? Yes Skin? Yes Ingestion? Yes
Routes of Entry:	Inhalation? Yes Skin? Yes Ingestion? Yes	Health Hazards (Acute and Chronic)	None
Health Hazards (Acute and Chronic)	None	Carcinogenicity:	None identified NTP? IARC Monographs? OSHA Regulation?
Carcinogenicity:	None identified NTP? IARC Monographs? OSHA Regulation?	Signs and Symptoms of Exposure	Irritation to upper respiratory tract, skin, eyes
Signs and Symptoms of Exposure	Irritation to upper respiratory tract, skin, eyes	Medical Conditions Generally Aggravated by Exposure	None
Medical Conditions Generally Aggravated by Exposure	None	Emergency First Aid Procedures	Ingestion: If conscious, give large amounts of water Eyes: Flush with water Inhalation: Move to fresh air Skin: Wash with soap and water
Emergency First Aid Procedures	Ingestion: If conscious, give large amounts of water Eyes: Flush with water Inhalation: Move to fresh air Skin: Wash with soap and water	Section VII - Precautions for Safe Handling and Use	Steps to be Taken in case Material is Released or Spilled: Wipe up spill and rinse with water, or collect in absorbent material and dispose of the absorbent material Waste Disposal Method Dispose in accordance with all applicable federal, state, and local environmental regulations. Precautions to be Taken in Handling and Storing Avoid eye and skin contact.
Section VII - Precautions for Safe Handling and Use	Steps to be Taken in case Material is Released or Spilled: Wipe up spill and rinse with water, or collect in absorbent material and dispose of the absorbent material Waste Disposal Method Dispose in accordance with all applicable federal, state, and local environmental regulations. Precautions to be Taken in Handling and Storing Avoid eye and skin contact.	Other Precautions	None
Other Precautions	None	Section VIII - Control Measures	Respiratory Protection (Specify Type) None
Section VIII - Control Measures	Respiratory Protection (Specify Type) None	Ventilation	Local Exhaust Yes Mechanical (General) Yes Special Other None
Ventilation	Local Exhaust Yes Mechanical (General) Yes Special Other None	Protective Gloves	Yes Eye Protection Safety goggles
Protective Gloves	Yes Eye Protection Safety goggles	Other Protective Clothing or Equipment	None
Other Protective Clothing or Equipment	None	Work/Hygiene Practices	None

Section V - Reactivity Data		Section V - Reactivity Data	
Stability	Unstable X Stable	Conditions to Avoid	None
Incompatibility	Strong oxidizing agents	Hazardous Decomposition or Byproducts	Carbon monoxide, Carbon dioxide
Hazardous Decomposition or Byproducts	Carbon monoxide, Carbon dioxide	Hazardous Polymerization	None
Hazardous Polymerization	None	Routes of Entry:	Inhalation? Yes Skin? Yes Ingestion? Yes
Routes of Entry:	Inhalation? Yes Skin? Yes Ingestion? Yes	Health Hazards (Acute and Chronic)	None
Health Hazards (Acute and Chronic)	None	Carcinogenicity:	None identified NTP? IARC Monographs? OSHA Regulation?
Carcinogenicity:	None identified NTP? IARC Monographs? OSHA Regulation?	Signs and Symptoms of Exposure	Irritation to upper respiratory tract, skin, eyes
Signs and Symptoms of Exposure	Irritation to upper respiratory tract, skin, eyes	Medical Conditions Generally Aggravated by Exposure	None
Medical Conditions Generally Aggravated by Exposure	None	Emergency First Aid Procedures	Ingestion: If conscious, give large amounts of water Eyes: Flush with water Inhalation: Move to fresh air Skin: Wash with soap and water
Emergency First Aid Procedures	Ingestion: If conscious, give large amounts of water Eyes: Flush with water Inhalation: Move to fresh air Skin: Wash with soap and water	Section VII - Precautions for Safe Handling and Use	Steps to be Taken in case Material is Released or Spilled: Wipe up spill and rinse with water, or collect in absorbent material and dispose of the absorbent material Waste Disposal Method Dispose in accordance with all applicable federal, state, and local environmental regulations. Precautions to be Taken in Handling and Storing Avoid eye and skin contact.
Section VII - Precautions for Safe Handling and Use	Steps to be Taken in case Material is Released or Spilled: Wipe up spill and rinse with water, or collect in absorbent material and dispose of the absorbent material Waste Disposal Method Dispose in accordance with all applicable federal, state, and local environmental regulations. Precautions to be Taken in Handling and Storing Avoid eye and skin contact.	Other Precautions	None
Other Precautions	None	Section VIII - Control Measures	Respiratory Protection (Specify Type) None
Section VIII - Control Measures	Respiratory Protection (Specify Type) None	Ventilation	Local Exhaust Yes Mechanical (General) Yes Special Other None
Ventilation	Local Exhaust Yes Mechanical (General) Yes Special Other None	Protective Gloves	Yes Eye Protection Safety goggles
Protective Gloves	Yes Eye Protection Safety goggles	Other Protective Clothing or Equipment	None
Other Protective Clothing or Equipment	None	Work/Hygiene Practices	None

Section V - Reactivity Data		Section V - Reactivity Data	
Stability	Unstable X Stable	Conditions to Avoid	None
Incompatibility	Copper, iron, silver salts, hydrogen peroxide, phenol, picric acid, formaldehyde, ether, alcohol, nitrogen oxide, strong bases, oxidizing agents.	Hazardous Decomposition or Byproducts	Carbon monoxide, carbon dioxide, nitrogen oxides, chlorine, hydrogen chloride
Hazardous Decomposition or Byproducts	Carbon monoxide, carbon dioxide, nitrogen oxides, chlorine, hydrogen chloride	Hazardous Polymerization	None
Hazardous Polymerization	None	Routes of Entry:	Inhalation? Yes Skin? Yes Ingestion? Yes
Routes of Entry:	Inhalation? Yes Skin? Yes Ingestion? Yes	Health Hazards (Acute and Chronic)	Toxicity has not been quantified. Sensitizing reactions (allergy) may occur by skin penetration including anaphylactic shock.
Health Hazards (Acute and Chronic)	Toxicity has not been quantified. Sensitizing reactions (allergy) may occur by skin penetration including anaphylactic shock.	Carcinogenicity:	None IARC Monographs? OSHA Regulation? No data
Carcinogenicity:	None IARC Monographs? OSHA Regulation? No data	Signs and Symptoms of Exposure	May cause irritation to skin, eyes, mucous membrane and upper respiratory tract.
Signs and Symptoms of Exposure	May cause irritation to skin, eyes, mucous membrane and upper respiratory tract.	Medical Conditions Generally Aggravated by Exposure	Respiratory conditions
Medical Conditions Generally Aggravated by Exposure	Respiratory conditions	Emergency First Aid Procedures	Treat symptomatically and supportively
Emergency First Aid Procedures	Treat symptomatically and supportively	Section VII - Precautions for Safe Handling and Use	Steps to be Taken in case Material is Released or Spilled: Wipe up spill with absorbent material. Dispose of properly. Waste Disposal Method Mix with vermiculite and dry caustic, wrap in paper and burn in a chemical incinerator equipped with afterburner and scrubber. Ignite in presence of sodium carbonate and slaked lime (CaOH) Precautions to be Taken in Handling and Storing Wear protective gear to avoid skin/eye contact.
Section VII - Precautions for Safe Handling and Use	Steps to be Taken in case Material is Released or Spilled: Wipe up spill with absorbent material. Dispose of properly. Waste Disposal Method Mix with vermiculite and dry caustic, wrap in paper and burn in a chemical incinerator equipped with afterburner and scrubber. Ignite in presence of sodium carbonate and slaked lime (CaOH) Precautions to be Taken in Handling and Storing Wear protective gear to avoid skin/eye contact.	Other Precautions	None
Other Precautions	None	Section VIII - Control Measures	Respiratory Protection (Specify Type) Chemical cartridge respirator with organic vapor cartridge.
Section VIII - Control Measures	Respiratory Protection (Specify Type) Chemical cartridge respirator with organic vapor cartridge.	Ventilation	Local Exhaust Yes Mechanical (General) No Special Other None
Ventilation	Local Exhaust Yes Mechanical (General) No Special Other None	Protective Gloves	Yes Eye Protection Splash proof goggles
Protective Gloves	Yes Eye Protection Splash proof goggles	Other Protective Clothing or Equipment	Protection to avoid skin contact
Other Protective Clothing or Equipment	Protection to avoid skin contact	Work/Hygiene Practices	Do not ingest. Avoid contact with skin, eyes and clothing. Wash thoroughly after handling.
Work/Hygiene Practices	Do not ingest. Avoid contact with skin, eyes and clothing. Wash thoroughly after handling.		

<p><b>EDVOTEK.</b></p> <p>Material Safety Data Sheet</p> <p>May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.</p>		<p>IDENTITY (As Used on Label and List)</p> <p>50% Glycerol</p>		<p>Notes: Blank spaces are not permitted. If any blank is not applicable, or no information is available, the space must be marked to indicate that.</p>	
<p><b>Section I</b></p> <p>Manufacturer's Name EDVOTEK, Inc.</p> <p>Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850</p> <p>Emergency Telephone Number (301) 251-5990 Telephone Number for information (301) 251-5990 Date Prepared 12/01/06 Signature of Preparer (optional)</p>		<p><b>Section II - Hazardous Ingredients/Identify Information</b></p> <p>Hazardous Components Specific Gravity (H<sub>2</sub>O = 1) % (Optional) Ethidium Bromide 0.0001% (0.0001%)</p> <p>This product contains no hazardous materials as defined by the OSHA Hazard Communication Standard.</p>		<p><b>Section III - Physical/Chemical Characteristics</b></p> <p>Boiling Point No data Specific Gravity (H<sub>2</sub>O = 1) No data</p> <p>Vapor Pressure (mm Hg.) No data Melting Point N/A</p> <p>Vapor Density (AIR = 1) No data Evaporation Rate (Boyl/Accurate = 1) No data</p> <p>Solubility in Water Soluble</p> <p>Appearance and Odor Blue liquid, no odor</p>	
<p><b>Section IV - Physical/Chemical Characteristics</b></p> <p>Flash Point (Method Used) No data</p> <p>Flammable Limits LEL No data UEL No data</p> <p>Extinguishing Media Dry chemical, carbon dioxide, water spray or foam</p> <p>Special Fire Fighting Procedures Use agents suitable for type of surrounding fire. Keep upwind, avoid breathing hazardous sulfur oxides and bromides. Wear SCBA.</p> <p>Unusual Fire and Explosion Hazards Unknown</p>		<p><b>Section V - Reactivity Data</b></p> <p>Stability Unstable X Stable None</p> <p>Incompatibility None known</p> <p>Hazardous Decomposition or Byproducts Sulfur oxides and bromides</p> <p>Hazardous Polymerization Will Not Occur X</p> <p><b>Section VI - Health Hazard Data</b></p> <p>Routes of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes</p> <p>Health Hazards (Acute and Chronic) May cause skin irritation, prolonged exposure can cause acute eye contact. May cause irritation. No data available for other routes.</p> <p>Carcinogenicity: NTP? No data IARC Monographs? OSHA Regulation? No data</p> <p>Signs and Symptoms of Exposure May cause skin or eye irritation</p> <p>Medical Conditions Generally Aggravated by Exposure None reported</p> <p>Emergency First Aid Procedures Treat symptomatically and supportively. Rinse contacted area with copious amounts of water.</p>		<p><b>Section VII - Control Measures</b></p> <p>Respiratory Protection (Specify Type) Chemical cartridge respirator with organic vapor cartridge.</p> <p>Ventilation Local Exhaust Yes Special None</p> <p>Protective Gloves Mechanical (General) Yes Other None</p> <p>Other Protective Clothing or Equipment Eye Protection Splash proof goggles</p> <p>Work/Hygiene Practices Do not ingest. Avoid contact with skin, eyes and clothing. Wash thoroughly after handling.</p>	

<p><b>EDVOTEK.</b></p> <p>Material Safety Data Sheet</p> <p>May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.</p>		<p>IDENTITY (As Used on Label and List)</p> <p>50% Glycerol</p>		<p>Notes: Blank spaces are not permitted. If any blank is not applicable, or no information is available, the space must be marked to indicate that.</p>	
<p><b>Section I</b></p> <p>Manufacturer's Name EDVOTEK, Inc.</p> <p>Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850</p> <p>Emergency Telephone Number (301) 251-5990 Telephone Number for information (301) 251-5990 Date Prepared 12/01/06 Signature of Preparer (optional)</p>		<p><b>Section II - Hazardous Ingredients/Identify Information</b></p> <p>Hazardous Components Specific Gravity (H<sub>2</sub>O = 1) % (Optional) CAS # 56-81-5</p>		<p><b>Section III - Physical/Chemical Characteristics</b></p> <p>Boiling Point 182°C Specific Gravity (H<sub>2</sub>O = 1) 1.26</p> <p>Vapor Pressure (mm Hg.) 20°C &lt;1MM Melting Point 20°C</p> <p>Vapor Density (AIR = 1) 3.1 Evaporation Rate (Boyl/Accurate = 1) N/A</p> <p>Solubility in Water Viscous colorless liquid</p> <p>Appearance and Odor Viscous colorless liquid</p>	
<p><b>Section IV - Physical/Chemical Characteristics</b></p> <p>Flash Point (Method Used) No data</p> <p>Flammable Limits LEL UEL</p> <p>Extinguishing Media Water spray, CO<sub>2</sub>, Dry chemical powder or appropriate foam</p> <p>Special Fire Fighting Procedures Wear SCBA and protective clothing to prevent contact with skin and eyes</p> <p>Unusual Fire and Explosion Hazards Emits toxic fumes under fire conditions.</p>		<p><b>Section V - Reactivity Data</b></p> <p>Stability Unstable X Stable None</p> <p>Incompatibility Strong oxidizing agents, strong bases, protect from moisture</p> <p>Hazardous Decomposition or Byproducts CO, CO<sub>2</sub></p> <p>Hazardous Polymerization Will Not Occur X</p> <p><b>Section VI - Health Hazard Data</b></p> <p>Routes of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes</p> <p>Health Hazards (Acute and Chronic) May cause skin irritation, prolonged exposure can cause headache, nausea and dizziness.</p> <p>Carcinogenicity: NTP? OSHA Regulation?</p> <p>Signs and Symptoms of Exposure Irritation to mucous membranes and upper respiratory tract</p> <p>Medical Conditions Generally Aggravated by Exposure None</p> <p>Emergency First Aid Procedures Flush with copious amounts of H<sub>2</sub>O for at least 15 minutes. Remove to fresh air &amp; remove contaminated clothing</p>		<p><b>Section VII - Control Measures</b></p> <p>Respiratory Protection (Specify Type) None</p> <p>Ventilation Local Exhaust X Mechanical (General) Mechanical exhaust Other</p> <p>Protective Gloves Chemical Resistant Eye Protection Safety goggles</p> <p>Other Protective Clothing or Equipment Wash thoroughly after handling</p>	

<p><b>EDVOTEK.</b></p> <p>Material Safety Data Sheet</p> <p>May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.</p>		<p>IDENTITY (As Used on Label and List)</p> <p>InstaStain® Ethidium Bromide</p>		<p>Notes: Blank spaces are not permitted. If any blank is not applicable, or no information is available, the space must be marked to indicate that.</p>	
<p><b>Section I</b></p> <p>Manufacturer's Name InstaStain, Inc.</p> <p>Address (Number, Street, City, State, Zip Code) P.O. Box 1232 West Bethesda, MD 20827</p> <p>Emergency Telephone Number (301) 251-5990 Telephone Number for information (301) 251-5990 Date Prepared 10/05/06 Signature of Preparer (optional)</p>		<p><b>Section II - Hazardous Ingredients/Identify Information</b></p> <p>Hazardous Components Specific Gravity (H<sub>2</sub>O = 1) % (Optional) Ethidium Bromide 0.0001% (0.0001%)</p> <p>(2,7-Diamino-10-Ethyl-9-Phenylphenanthridinium Bromide)</p> <p>CAS# 139-33-3</p>		<p><b>Section III - Physical/Chemical Characteristics</b></p> <p>Boiling Point No data Specific Gravity (H<sub>2</sub>O = 1) No data</p> <p>Vapor Pressure (mm Hg.) No data Melting Point No data</p> <p>Vapor Density (AIR = 1) No data Evaporation Rate (Boyl/Accurate = 1) No data</p> <p>Solubility in Water Soluble</p> <p>Appearance and Odor Chemical bound to paper, no odor</p>	
<p><b>Section IV - Physical/Chemical Characteristics</b></p> <p>Flash Point (Method Used) No data</p> <p>Flammable Limits LEL N.D. UEL N.D.</p> <p>Extinguishing Media Water spray, carbon dioxide, dry chemical powder, alcohol or polymer foam</p> <p>Special Fire Fighting Procedures Wear protective clothing and SCBA to prevent contact with skin &amp; eyes</p> <p>Unusual Fire and Explosion Hazards Emits toxic fumes</p>		<p><b>Section V - Reactivity Data</b></p> <p>Stability Unstable X Stable None</p> <p>Incompatibility Strong oxidizing agents</p> <p>Hazardous Decomposition or Byproducts Carbon dioxide, nitrogen oxide, hydrogen bromide gas</p> <p>Hazardous Polymerization Will Not Occur X</p> <p><b>Section VI - Health Hazard Data</b></p> <p>Routes of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes</p> <p>Health Hazards (Acute and Chronic) May alter genetic material. Acute: Material irritating to mucous membranes, upper respiratory tract, eyes, skin</p> <p>Carcinogenicity: No data available NTP? IARC Monographs? OSHA Regulation?</p> <p>Signs and Symptoms of Exposure Irritation to mucous membranes and upper respiratory tract</p> <p>Medical Conditions Generally Aggravated by Exposure No data</p> <p>Emergency First Aid Procedures Treat symptomatically and supportively</p>		<p><b>Section VII - Control Measures</b></p> <p>Respiratory Protection (Specify Type) SCBA</p> <p>Ventilation Local Exhaust Yes Special Chem. fume hood</p> <p>Protective Gloves Rubber Eye Protection Chem. safety goggles</p> <p>Other Protective Clothing or Equipment Rubber boots</p> <p>Work/Hygiene Practices Use in chemical fume hood with proper protective lab gear.</p>	

<p><b>EDVOTEK.</b></p> <p>Material Safety Data Sheet</p> <p>May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.</p>		<p>IDENTITY (As Used on Label and List)</p> <p>InstaStain® Ethidium Bromide</p>		<p>Notes: Blank spaces are not permitted. If any blank is not applicable, or no information is available, the space must be marked to indicate that.</p>	
<p><b>Section I</b></p> <p>Manufacturer's Name InstaStain, Inc.</p> <p>Address (Number, Street, City, State, Zip Code) P.O. Box 1232 West Bethesda, MD 20827</p> <p>Emergency Telephone Number (301) 251-5990 Telephone Number for information (301) 251-5990 Date Prepared 10/05/06 Signature of Preparer (optional)</p>		<p><b>Section II - Hazardous Ingredients/Identify Information</b></p> <p>Hazardous Components Specific Gravity (H<sub>2</sub>O = 1) % (Optional) Ethidium Bromide 0.0001% (0.0001%)</p> <p>(2,7-Diamino-10-Ethyl-9-Phenylphenanthridinium Bromide)</p> <p>CAS# 139-33-3</p>		<p><b>Section III - Physical/Chemical Characteristics</b></p> <p>Boiling Point No data Specific Gravity (H<sub>2</sub>O = 1) No data</p> <p>Vapor Pressure (mm Hg.) No data Melting Point No data</p> <p>Vapor Density (AIR = 1) No data Evaporation Rate (Boyl/Accurate = 1) No data</p> <p>Solubility in Water Soluble</p> <p>Appearance and Odor Chemical bound to paper, no odor</p>	
<p><b>Section IV - Physical/Chemical Characteristics</b></p> <p>Flash Point (Method Used) No data</p> <p>Flammable Limits LEL N.D. UEL N.D.</p> <p>Extinguishing Media Water spray, carbon dioxide, dry chemical powder, alcohol or polymer foam</p> <p>Special Fire Fighting Procedures Wear protective clothing and SCBA to prevent contact with skin &amp; eyes</p> <p>Unusual Fire and Explosion Hazards Emits toxic fumes</p>		<p><b>Section V - Reactivity Data</b></p> <p>Stability Unstable X Stable None</p> <p>Incompatibility Strong oxidizing agents</p> <p>Hazardous Decomposition or Byproducts Carbon dioxide, nitrogen oxide, hydrogen bromide gas</p> <p>Hazardous Polymerization Will Not Occur X</p> <p><b>Section VI - Health Hazard Data</b></p> <p>Routes of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes</p> <p>Health Hazards (Acute and Chronic) May alter genetic material. Acute: Material irritating to mucous membranes, upper respiratory tract, eyes, skin</p> <p>Carcinogenicity: No data available NTP? IARC Monographs? OSHA Regulation?</p> <p>Signs and Symptoms of Exposure Irritation to mucous membranes and upper respiratory tract</p> <p>Medical Conditions Generally Aggravated by Exposure No data</p> <p>Emergency First Aid Procedures Treat symptomatically and supportively</p>		<p><b>Section VII - Control Measures</b></p> <p>Respiratory Protection (Specify Type) SCBA</p> <p>Ventilation Local Exhaust Yes Special Chem. fume hood</p> <p>Protective Gloves Rubber Eye Protection Chem. safety goggles</p> <p>Other Protective Clothing or Equipment Rubber boots</p> <p>Work/Hygiene Practices Use in chemical fume hood with proper protective lab gear.</p>	

<p><b>EDVOTEK.</b></p> <p>Material Safety Data Sheet</p> <p>May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.</p>		<p>IDENTITY (As Used on Label and List)</p> <p>InstaStain® Ethidium Bromide</p>		<p>Notes: Blank spaces are not permitted. If any blank is not applicable, or no information is available, the space must be marked to indicate that.</p>	
<p><b>Section I</b></p> <p>Manufacturer's Name InstaStain, Inc.</p> <p>Address (Number, Street, City, State, Zip Code) P.O. Box 1232 West Bethesda, MD 20827</p> <p>Emergency Telephone Number (301) 251-5990 Telephone Number for information (301) 251-5990 Date Prepared 10/05/06 Signature of Preparer (optional)</p>		<p><b>Section II - Hazardous Ingredients/Identify Information</b></p> <p>Hazardous Components Specific Gravity (H<sub>2</sub>O = 1) % (Optional) Ethidium Bromide 0.0001% (0.0001%)</p> <p>(2,7-Diamino-10-Ethyl-9-Phenylphenanthridinium Bromide)</p> <p>CAS# 139-33-3</p>		<p><b>Section III - Physical/Chemical Characteristics</b></p> <p>Boiling Point No data Specific Gravity (H<sub>2</sub>O = 1) No data</p> <p>Vapor Pressure (mm Hg.) No data Melting Point No data</p> <p>Vapor Density (AIR = 1) No data Evaporation Rate (Boyl/Accurate = 1) No data</p> <p>Solubility in Water Soluble</p> <p>Appearance and Odor Chemical bound to paper, no odor</p>	
<p><b>Section IV - Physical/Chemical Characteristics</b></p> <p>Flash Point (Method Used) No data</p> <p>Flammable Limits LEL N.D. UEL N.D.</p> <p>Extinguishing Media Water spray, carbon dioxide, dry chemical powder, alcohol or polymer foam</p> <p>Special Fire Fighting Procedures Wear protective clothing and SCBA to prevent contact with skin &amp; eyes</p> <p>Unusual Fire and Explosion Hazards Emits toxic fumes</p>		<p><b>Section V - Reactivity Data</b></p> <p>Stability Unstable X Stable None</p> <p>Incompatibility Strong oxidizing agents</p> <p>Hazardous Decomposition or Byproducts Carbon dioxide, nitrogen oxide, hydrogen bromide gas</p> <p>Hazardous Polymerization Will Not Occur X</p> <p><b>Section VI - Health Hazard Data</b></p> <p>Routes of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes</p> <p>Health Hazards (Acute and Chronic) May alter genetic material. Acute: Material irritating to mucous membranes, upper respiratory tract, eyes, skin</p> <p>Carcinogenicity: No data available NTP? IARC Monographs? OSHA Regulation?</p> <p>Signs and Symptoms of Exposure Irritation to mucous membranes and upper respiratory tract</p> <p>Medical Conditions Generally Aggravated by Exposure No data</p> <p>Emergency First Aid Procedures Treat symptomatically and supportively</p>		<p><b>Section VII - Control Measures</b></p> <p>Respiratory Protection (Specify Type) SCBA</p> <p>Ventilation Local Exhaust Yes Special Chem. fume hood</p> <p>Protective Gloves Rubber Eye Protection Chem. safety goggles</p> <p>Other Protective Clothing or Equipment Rubber boots</p> <p>Work/Hygiene Practices Use in chemical fume hood with proper protective lab gear.</p>	

<p><b>EDVOTEK.</b></p> <p>Material Safety Data Sheet</p> <p>May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.</p>		<p>IDENTITY (As Used on Label and List)</p> <p>InstaStain® Ethidium Bromide</p>		<p>Notes: Blank spaces are not permitted. If any blank is not applicable, or no information is available, the space must be marked to indicate that.</p>	
<p><b>Section I</b></p> <p>Manufacturer's Name InstaStain, Inc.</p> <p>Address (Number, Street, City, State, Zip Code) P.O. Box 1232 West Bethesda, MD 20827</p> <p>Emergency Telephone Number (301) 251-5990 Telephone Number for information (301) 251-5990 Date Prepared 10/05/06 Signature of Preparer (optional)</p>		<p><b>Section II - Hazardous Ingredients/Identify Information</b></p> <p>Hazardous Components Specific Gravity (H<sub>2</sub>O = 1) % (Optional) Ethidium Bromide 0.0001% (0.0001%)</p> <p>(2,7-Diamino-10-Ethyl-9-Phenylphenanthridinium Bromide)</p> <p>CAS# 139-33-3</p>		<p><b>Section III - Physical/Chemical Characteristics</b></p> <p>Boiling Point No data Specific Gravity (H<sub>2</sub>O = 1) No data</p> <p>Vapor Pressure (mm Hg.) No data Melting Point No data</p> <p>Vapor Density (AIR = 1) No data Evaporation Rate (Boyl/Accurate = 1) No data</p> <p>Solubility in Water Soluble</p> <p>Appearance and Odor Chemical bound to paper, no odor</p>	
<p><b>Section IV - Physical/Chemical Characteristics</b></p> <p>Flash Point (Method Used) No data</p> <p>Flammable Limits LEL N.D. UEL N.D.</p> <p>Extinguishing Media Water spray, carbon dioxide, dry chemical powder, alcohol or polymer foam</p> <p>Special Fire Fighting Procedures Wear protective clothing and SCBA to prevent contact with skin &amp; eyes</p> <p>Unusual Fire and Explosion Hazards Emits toxic fumes</p>		<p><b>Section V - Reactivity Data</b></p> <p>Stability Unstable X Stable None</p> <p>Incompatibility Strong oxidizing agents</p> <p>Hazardous Decomposition or Byproducts Carbon dioxide, nitrogen oxide, hydrogen bromide gas</p> <p>Hazardous Polymerization Will Not Occur X</p> <p><b>Section VI - Health Hazard Data</b></p> <p>Routes of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes</p> <p>Health Hazards (Acute and Chronic) May alter genetic material. Acute: Material irritating to mucous membranes, upper respiratory tract, eyes, skin</p> <p>Carcinogenicity: No data available NTP? IARC Monographs? OSHA Regulation?</p> <p>Signs and Symptoms of Exposure Irritation to mucous membranes and upper respiratory tract</p> <p>Medical Conditions Generally Aggravated by Exposure No data</p> <p>Emergency First Aid Procedures Treat symptomatically and supportively</p>		<p><b>Section VII - Control Measures</b></p> <p>Respiratory Protection (Specify Type) SCBA</p> <p>Ventilation Local Exhaust Yes Special Chem. fume hood</p> <p>Protective Gloves Rubber Eye Protection Chem. safety goggles</p> <p>Other Protective Clothing or Equipment Rubber boots</p> <p>Work/Hygiene Practices Use in chemical fume hood with proper protective lab gear.</p>	



**Material Safety Data Sheet**

May be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.

**IDENTITY** (As Used on Label and LHD)  
Potassium Chloride

**Section I**

Manufacturer's Name  
EDVOTEK, Inc.  
Address (Number, Street, City, State, Zip Code)  
14676 Rothgeb Drive  
Rockville, MD 20850

Emergency Telephone Number (301) 251-5990  
Telephone Number for information (301) 251-5990  
Date Prepared 1/20/06  
Signature of Preparer (optional)

Note: Blank space not identified. If any item is not applicable, or no information is available, the space must be marked to indicate that.

**Section II - Hazardous Ingredients/Identify Information**

Hazardous Components (Specific Chemical Identity, Common Name(s))  
OSHA PEL: N/A  
ACGIH TLV: N/A  
Recombined % (Optional): N/A

**Section III - Physical/Chemical Characteristics**

Boiling Point: 1500°C  
Specific Gravity (H<sub>2</sub>O = 1): 1.984  
Vapor Pressure (mm Hg): N/A  
Melting Point: N/A  
Vapor Density (AIR = 1): N/A  
Evaporation Rate (Easy/Accurate = 1): N/A  
Solubility in Water: Soluble  
Appearance and Odor: White crystalline powder, odorless

**Section IV - Physical/Chemical Characteristics**

Flash Point (Method Used): N/A  
Flammable Limits: LEL: N/A, UEL: N/A  
Extinguishing Media: Dry chemical, carbon dioxide, water spray or foam.  
Special Fire Fighting Procedures: Wear SCBA and protective clothing  
Unusual Fire and Explosion Hazards: None known

**Section V - Reactivity Data**

Stability: Unstable: X, Stable: X  
Conditions to Avoid: Violent reaction w/ BrF<sub>3</sub>, (H<sub>2</sub>SO<sub>4</sub>+KMNO<sub>4</sub>)  
Incompatibility: BrF<sub>3</sub>

Hazardous Decomposition or Byproducts: N/A  
Hazardous Polymerization: May Occur: X, Will Not Occur: X  
Conditions to Avoid: N/A

**Section VI - Health Hazard Data**

Route(s) of Entry: Inhalation? Yes, Skin? Yes, Ingestion? Yes  
Health Hazards (Acute and Chronic): Presents no serious health effects in normal industrial uses. However, dust may cause irritation of eyes.  
Carcinogenicity: NTP? No data, IARC? No data, OSHA Regulation?  
Signs and Symptoms of Exposure: Irritation of eyes, skin and mucous membranes, vomiting, weakness and circulatory problems.  
Medical Conditions Generally Aggravated by Exposure: Unknown

**Section VII - Precautions for Safe Handling and Use**

Emergency First Aid Procedures: Eyes and skin: Flush with copious amounts of water for 15 min. Inhalation: move to fresh air. Ingestion: get medical attention. Never give anything by mouth to an unconscious person.  
Steps to be Taken in case Material is Released or Spilled: Sweep up material and place in to suitable waste container. Flush spill area with water.  
Waste Disposal Method: Dispose of material in accordance with state, local and Federal regulations.  
Precautions to be Taken in Handling and Storage: Store in tightly sealed container. Protect from moisture.  
Other Precautions: Avoid skin contact and avoid breathing dust.

**Section VIII - Control Measures**

Respiratory Protection (Specify Type): NIOSH/MSHA approved respirator

Ventilation	Local Exhaust	Recommended	Special
	Mechanical (General)	Yes	Other
			N/A

Protective Gloves: PE or rubber gloves  
Eye Protection: Safety glasses or goggles  
Other Protective Clothing or Equipment: Lab coat  
Work/Hygiene Practices: Wash thoroughly after use.

**Material Safety Data Sheets**

Full size (8.5 x 11") pdf copy of MSDS available at [www.edvotek.com](http://www.edvotek.com) or by request.

EDVO-Kit #  
**302**